LANAFORM®

Office Massage Chair

Massager seat

Your perfect partner

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Uw beste bondgenoot

OFFICE MASSAGE CHAIR

USER MANUAL

Thank you for choosing the OFFICE MASSAGE CHAIR from Lanaform®, a high-quality massage device. The OFFICE MASSAGE CHAIR is a massaging office chair equipped with four massage heads which help reduce tension and muscular pain and soothe tired muscles. Comfortable and sophisticated, the OFFICE MASSAGE CHAIR offers a very pleasant 'shiatsu' and 'rolling' massage all over your back. This appliance is also equipped with a vibration function for the seat and an automatic function offering a pleasant blend of the different massage options. What's more, it is operated by an easy-to-use control which is directly integrated into the armrest of your chair. The OFFICE MASSAGE CHAIR can be used in the office, at home, etc.

PLEASE READ ALL THE INSTRUCTIONS PRIOR TO USING YOUR OFFICE MASSAGE CHAIR, IN PARTICULAR THESE FEW BASIC SAFETY INSTRUCTIONS:

- Only use this appliance for the purposes described in this manual.
- Before using the device, remember to remove the silica sachet used to combat humidity in the air thus protecting the packaged items.
- This appliance is not intended for use by persons (including children) whose physical, sensory or mental capabilities are impaired, or by persons lacking experience or knowledge, unless they have had the benefit of prior supervision and instruction in using the device from a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
- Please do not use accessories which are not recommended by LANAFORM® or which are not supplied with this appliance.
- If the supply cord is damaged, it must be replaced by a special cord or assembly available from the manufacturer or its service agent.
- Do not use this appliance if the plug is damaged, if it is not working correctly, if it has been dropped or if it is damaged or if it has fallen into water. If this is the case, have the appliance inspected and repaired by the supplier or its after-sales department.
- Do not carry this appliance by its supply lead or use this lead as a handle.
- Never use drawing pins or other metal attachments with this appliance.
- Always disconnect the appliance from the mains after use or before cleaning it.
- Should the appliance fall into water, unplug it immediately before retrieving it.
- An electrical appliance should never be left plugged in and unsupervised. Disconnect it from the mains if you are not using it.
- Keep the power lead away from hot surfaces.
- Do not use this product before going to bed. Massage has a stimulating effect which may inhibit sleep.
- Do not use this massage appliance while you are asleep.
- Never let anything fall into and never insert anything into one of the openings.
- Never use this appliance in a room where aerosol products (sprays) are used or in a room where oxygen is administered.
- Do not use this appliance under a cover or under a cushion. Excessive heat may cause a fire, electrocution of the individual or other injuries.
- When unplugging the appliance, please make sure that the  button is switched off before removing the plug.
- If you are worried about your health, consult a doctor before using this appliance.
- If you feel any pain while using this appliance, stop using it immediately and consult your doctor.
- The appliance has a heated surface. Person insensitive to heat must be careful when using the appliance.
- It is not recommended to use the device for more than 15 minutes at a time.
- Never use this appliance on a part of the body which is swollen or inflamed or if you have any form of skin rash.
- This product is a non-professional massage appliance designed to soothe tired muscles. Do not use this appliance as a replacement for medical care, it is intended exclusively for family use.

HOW TO USE

1) Details :

- 4 massage heads
- 'shiatsu' & 'rolling' massage
- Adjustable massage head spacing
- 1 'auto' mode
- 3 massage areas (entire back, upper back and lower back)
- 1 "targeted massage" function
- Vibrating seat
- "Heating" function
- Retractable flex on the side of the chair

2) Description of the OFFICE MESSAGE CHAIR :

Control device:

- | | |
|----------------------|-------------------------|
| ① On/Off | ⑦ Massage head spacing* |
| ② "Auto" mode | ⑧ Full back massage |
| ③ "Shiatsu" function | ⑨ Upper back massage |
| ④ "Rolling" function | ⑩ Lower back massage |
| ⑤ "Heating" function | ⑪ Targeted massage |
| ⑥ Vibrating seat | |

* only for the "rolling" function.

USING OF THE OFFICE MESSAGE CHAIR

The OFFICE MESSAGE CHAIR has a retractable flex located on the right site of the chair.

Plug this into a power socket. Once connected to the power socket, the OFFICE MESSAGE CHAIR is ready for use.

On/Off :

Starting or stopping the appliance is controlled by simply pressing the red button (diagram, point 1).

Back:

Automatic mode (diagram, point 2):

This mode consists of a sequence of different massages. This function is valuable since it offers a varied massage.


Manual mode:

- You have the option of choosing between two types of massage: 'shiatsu' (diagram, point 3) and 'rolling' (diagram, point 4) for which you may choose to adjust the massage head spacing (diagram, point 7).
- You also have the option to select the location of your massage:
 - a) Targeted massage, precise location (diagram, point 11).
 - b) Your lower back (diagram, point 10).
 - c) Your upper back (diagram, point 9).
 - d) Entire back (diagram, point 8).
- The 'heating' function can be switched on at any time (diagram, point 5).

Warning: You may not feel the heat through clothing such as a jumper or jacket.

Seat:

- You also have the benefit of a vibrating seat that you can choose to enable or not (diagram, point 6).

MAINTENANCE

Storage

- Put the OFFICE MESSAGE CHAIR in a dry place.
- Avoid it coming into contact with sharp edges and other pointed objects which might cut or pierce the surface.
- To prevent the flex from breaking, do not move the armchair by pulling on it.

Cleaning

- Unplug the massaging office armchair at the power socket and wait a few minutes before cleaning it.
- Clean it with a soft, slightly damp sponge.
- Never let water or any other liquid come into contact with the chair.
- Never use abrasive detergents, brushes or a thinning agent to clean it.

WASTE MANAGEMENT

- All the packaging is composed of materials that pose no hazard for the environment and which can be disposed of at your local sorting centre to be used as secondary raw materials. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.
- When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law.

LIMITED WARRANTY

LANAFORM® guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, except in the circumstances described below.

The LANAFORM® guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this LANAFORM® product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorized accessories, changes made to the product or any other circumstance, of whatever sort, that is outside LANAFORM®'s control.

LANAFORM® may not be held liable for any type of circumstantial, indirect or specific damage.

All implicit guarantees relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof of purchase can be supplied.

On receipt, LANAFORM® will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM®'s Service Centre®. Any attempt to maintain this product by a person other than LANAFORM®'s Service Centre will render this guarantee void.

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